



ANSABERE MULTI-PURPOSE RETRACTOR Mod. 8020



Urology Kit Mod. 8020-U2



UNIT	REF.	DESCRIPTION	
A	1	B150	Rail clamp
B	1	B400	Toothed clamp
C	1	R300	Frame support bar
D	1	F300	Dopple socket joint frame
E	1	F150	Extension bridge
F	6	B502	Universal clamp
G	4	P205	Blade holder

UNIT	REF.	DESCRIPTION	
H	2	VC312M	Wall blade with handle 90x60mm
I	1	VP452	Malleable blade 50x170mm
J	1	VP462	Malleable blade 80x170mm
K	1	VP124	Narrow flat blade 40x180mm
L	1	VP258	Wide flat blade (Fence) 100x170mm
M	1	VP502	"V" shaped blade 72x130mm