




ANSABERE MULTI-PURPOSE RETRACTOR Mod. 8020



Vascular Surgery Kit Mod. 8020-V3



UNIT	REF.	DESCRIPTION
A	1	B150 Rail clamp
B	1	B400 Toothed clamp
C	1	R300 Frame support bar
D	1	F300 Doppie socket joint frame
E	1	F150 Extension bridge
F	8	B502 Universal clamp
G	3	P205 Blade holder
H	1	P405 Angling blade holder
I	1	O105 Flexible illumination terminal

UNIT	REF.	DESCRIPTION
J	2	VC312M Wall blade with handle 90x60mm
K	1	VC316M Wall blade with handle 100x60mm
L	1	VP452 Malleable blade 50x170mm
M	1	VP462 Malleable blade 80x170mm
N	1	VP140 Narrow flat blade 50x150mm
O	1	VP204 Medium flat blade 65x170mm
P	1	VP258 Wide flat blade (Fence) 100x170mm
Q	1	VG202 Renal vein blade 25x180mm