




RETRACTOR UNIVERSAL ANSABERE - Mod. 8020



Kit Cirugía Vasculr Mod. 8020-V3



UNID.	REF.	DESCRIPCIÓN	UNID.	REF.	DESCRIPCIÓN		
A	1	B150	Brida vertical	J	2	VC312M	Valva de pared con mango 90x60mm
B	1	B400	Brida dentada	K	1	VC316M	Valva de pared con mango 100x60mm
C	1	R300	Barra porta marco	L	1	VP452	Valva ranurada maleable 50x170mm
D	1	F300	Marco doble rótula	M	1	VP462	Valva ranurada maleable 80x170mm
E	1	F150	Accesorio de cierre	N	1	VP140	Valva plana estrecha 50x150mm
F	8	B502	Brida universal	O	1	VP204	Valva plana media 65x170mm
G	3	P205	Porta valvas	P	1	VP258	Valva plana ancha (Fence) 100x170mm
H	1	P405	Porta valvas basculante	Q	1	VG202	Valva renal 25x180mm
I	1	O105	Terminal iluminación flexible				