



ANSABERE MULTI-PURPOSE RETRACTOR



MULTI-PURPOSE RETRACTOR - Mod. 8020


- Independent frame arms
- Handle to fasten each frame arms

UNIT	REF.	DESCRIPTION
1	B150	Rail clamp
1	B400	Toothed clamp
1	R300	Frame support bar
1	F300	Dopple socket joint frame





EVOLUTION RETRACTOR SYSTEM - Mod. 8030

UNIT	REF.	DESCRIPTION
1	B152	Rail clamp
1	B402	Toothed clamp with double handle
1	R302	Frame support bar
1	F320	Detachable frame
2	F325	Frame arm with socket joint




B502 Universal clamp



P205 Blade holder

Blades assembling detail:




P405 Angling blade holder

Movement detail:



Allows the blade to be angled by rotating the handle



O105 Flexible illumination terminal

Includes connectors to attach to any light cable.
Autoclavable (134 °C / 273 °F)



Assembling detail:

F150 Extension bridge



Allows to close or complete the frame.

BLADES

Wall blade with handle

VC300-M	55x40 mm
VC308-M	70x50 mm
VC312-M	90x60 mm
VC316-M	100x60 mm
VC322-M	110x70 mm
VC324-M	110x80 mm



Wall blade

VC300	55x40 mm
VC308	70x50 mm
VC312	90x60 mm
VC316	100x60 mm
VC322	110x70 mm
VC324	110x80 mm



Narrow flat blade

VP084	25x130 mm
VP100	30x80 mm
VP104	30x120 mm
VP105	30x160 mm
VP122	40x130 mm
VP124	40x180 mm
VP140	50x150 mm
VP144	50x190 mm



Medium flat blade

VP200	65x110 mm
VP202	65x140 mm
VP204	65x170 mm
VP210	80x140 mm
VP212	80x160 mm
VP214	80x190 mm



Wide flat blade (Fence)

VP250	100x110 mm
VP254	100x130 mm
VP256	100x150 mm
VP258	100x170 mm
VP262	120x190 mm



Malleable blade

VP450	50x150 mm
VP452	50x170 mm
VP454	50x190 mm
VP460	80x150 mm
VP462	80x170 mm
VP464	80x190 mm



Pelvic blade (St. Mark's)

VP350	56x170 mm
VP352	56x190 mm
VP370	70x180 mm
VP372	70x200 mm



"V" - shaped blade

VP500	72x110 mm
VP502	72x130 mm
VP512	82x130 mm
VP514	82x150 mm



Harrington blade

VP600	50x130 mm
VP610	65x150 mm
VP614	65x180 mm
VP620	80x200 mm



Renal vein blade

VG192	25x80 mm
VG194	25x100 mm
VG196	25x120 mm
VG198	25x140 mm
VG200	25x150 mm
VG202	25x180 mm



Kelly blade

VP800	48x80 mm
VP802	48x100 mm
VP804	48x120 mm
VP810	60x150 mm
VP812	60x180 mm



Doyen blade

VP850	48x80 mm
VP852	48x100 mm
VP854	48x120 mm
VP860	60x150 mm
VP862	60x180 mm



Deaver blade

VC100	36x50 mm
VC120	60x60 mm
VC130	80x70 mm



Rectangular flat blade

VP300	85x60 mm
VP302	85x80 mm
VP304	85x100 mm



Fenestrated blade

VP700	28x60 mm
VP702	28x80 mm
VP704	28x100 mm
VP710	48x80 mm
VP712	48x100 mm
VP714	48x120 mm



Two fingers malleable blade

VS200	40 mm (width)
VS210	60 mm



Three fingers malleable blade

VS300	90 mm
VS310	120 mm



Four fingers malleable blade

VS400	120 mm
VS410	150 mm



Blades with handle

All blades are available with rotating handle adding an "M" to the reference

Example: VP124 <> VP124-M

