



# ANSABERE MULTI-PURPOSE RETRACTOR



## MULTI-PURPOSE RETRACTOR - Mod. 8020


- Independent frame arms
- Handle to fasten each frame arms

UNIT	REF.	DESCRIPTION
1	B150	Rail clamp
1	B400	Toothed clamp
1	R300	Frame support bar
1	F300	Dopple socket joint frame





## EVOLUTION RETRACTOR SYSTEM - Mod. 8030

UNIT	REF.	DESCRIPTION
1	B152	Rail clamp
1	B402	Toothed clamp with double handle
1	R302	Frame support bar
1	F320	Detachable frame
2	F325	Frame arm with socket joint




**B502 Universal clamp**



**P205 Blade holder**

Blades assembling detail:




**P405 Angling blade holder**

Movement detail:



Allows the blade to be angled by rotating the handle



**O105 Flexible illumination terminal**

Includes connectors to attach to any light cable.  
Autoclavable (134 °C / 273 °F)



Assembling detail:

**F150 Extension bridge**



Allows to close or complete the frame.

# BLADES

## Wall blade with handle

<b>VC300-M</b>	55x40 mm
<b>VC308-M</b>	70x50 mm
<b>VC312-M</b>	90x60 mm
<b>VC316-M</b>	100x60 mm
<b>VC322-M</b>	110x70 mm
<b>VC324-M</b>	110x80 mm



## Wall blade

<b>VC300</b>	55x40 mm
<b>VC308</b>	70x50 mm
<b>VC312</b>	90x60 mm
<b>VC316</b>	100x60 mm
<b>VC322</b>	110x70 mm
<b>VC324</b>	110x80 mm



## Narrow flat blade

<b>VP084</b>	25x130 mm
<b>VP100</b>	30x80 mm
<b>VP104</b>	30x120 mm
<b>VP105</b>	30x160 mm
<b>VP122</b>	40x130 mm
<b>VP124</b>	40x180 mm
<b>VP140</b>	50x150 mm
<b>VP144</b>	50x190 mm



## Medium flat blade

<b>VP200</b>	65x110 mm
<b>VP202</b>	65x140 mm
<b>VP204</b>	65x170 mm
<b>VP210</b>	80x140 mm
<b>VP212</b>	80x160 mm
<b>VP214</b>	80x190 mm



## Wide flat blade (Fence)

<b>VP250</b>	100x110 mm
<b>VP254</b>	100x130 mm
<b>VP256</b>	100x150 mm
<b>VP258</b>	100x170 mm
<b>VP262</b>	120x190 mm



## Malleable blade

<b>VP450</b>	50x150 mm
<b>VP452</b>	50x170 mm
<b>VP454</b>	50x190 mm
<b>VP460</b>	80x150 mm
<b>VP462</b>	80x170 mm
<b>VP464</b>	80x190 mm



## Pelvic blade (St. Mark's)

<b>VP350</b>	56x170 mm
<b>VP352</b>	56x190 mm
<b>VP370</b>	70x180 mm
<b>VP372</b>	70x200 mm



## "V" - shaped blade

<b>VP500</b>	72x110 mm
<b>VP502</b>	72x130 mm
<b>VP512</b>	82x130 mm
<b>VP514</b>	82x150 mm



## Harrington blade

<b>VP600</b>	50x130 mm
<b>VP610</b>	65x150 mm
<b>VP614</b>	65x180 mm
<b>VP620</b>	80x200 mm



## Renal vein blade

<b>VG192</b>	25x80 mm
<b>VG194</b>	25x100 mm
<b>VG196</b>	25x120 mm
<b>VG198</b>	25x140 mm
<b>VG200</b>	25x150 mm
<b>VG202</b>	25x180 mm



### Kelly blade

<b>VP800</b>	48x80 mm
<b>VP802</b>	48x100 mm
<b>VP804</b>	48x120 mm
<b>VP810</b>	60x150 mm
<b>VP812</b>	60x180 mm



### Doyen blade

<b>VP850</b>	48x80 mm
<b>VP852</b>	48x100 mm
<b>VP854</b>	48x120 mm
<b>VP860</b>	60x150 mm
<b>VP862</b>	60x180 mm



### Deaver blade

<b>VC100</b>	36x50 mm
<b>VC120</b>	60x60 mm
<b>VC130</b>	80x70 mm



### Rectangular flat blade

<b>VP300</b>	85x60 mm
<b>VP302</b>	85x80 mm
<b>VP304</b>	85x100 mm



### Fenestrated blade

<b>VP700</b>	28x60 mm
<b>VP702</b>	28x80 mm
<b>VP704</b>	28x100 mm
<b>VP710</b>	48x80 mm
<b>VP712</b>	48x100 mm
<b>VP714</b>	48x120 mm



### Two fingers malleable blade

<b>VS200</b>	40 mm (width)
<b>VS210</b>	60 mm



### Three fingers malleable blade

<b>VS300</b>	90 mm
<b>VS310</b>	120 mm



### Four fingers malleable blade

<b>VS400</b>	120 mm
<b>VS410</b>	150 mm



### Blades with handle

All blades are available with rotating handle adding an "M" to the reference

Example: **VP124 <> VP124-M**

